

**RIH - WRIST/HAND CT  
GE LIGHTSPEED 16 / OPTIMA CT580 PROTOCOL**

**Indication: fracture, dislocation, osteomyelitis, bone injury, bone tumor.**

<b>Position/Landmark</b>	Supine , feet first Zero Appropriately
<b>Topogram Direction</b>	Craniocaudal
<b>Respiratory Phase</b>	Any
<b>Scan Type</b>	Helical
<b>KV / mA / Rotation time (sec) Pitch / Speed (mm/rotation) Noise Index</b>	120kv / smart mA (100-440) / .5 sec .938:1 , 9.37mm 25.00
<b>Detector width x Rows = Beam Collimation</b>	0.625mm x 16 = 10mm
<b>Helical Set</b>	
Slice Thickness/ Spacing	recon      body      thickness/ part      spacing      recon
Algorithm	1   thin wrist/hand   .6mm x .6mm      bone      for dmpr
Recon Destination	2 <b>wrist/hand bone</b> 1.25mm x 1.25mm      bone      pacs
	3 <b>wrist/hand soft tissue</b> 1.25mm x 1.25 mm      standard      pacs
<b>Scan Start / End Locations</b>	determined by technologist or radiologist to include the anatomy of interest
<b>DFOV</b>	18cm decrease appropriately
<b>IV Contrast Volume / Type / Rate</b>	70cc omni 350 / 2cc per second if needed
<b>Scan Delay</b>	65 seconds
<b>Archiving to MOD</b>	Only prospective recons will be archived to mod as done by the scanner.
<b>2D/3D Technique Used</b>	DMPR of 2mm x 2mm coronal and sagittal wrist/hand series (auto-batch off), average mode, auto-transferred to PACS Also, there is a 2mm x 2mm true axial reformat if needed due to the patient's position. For wrist ct with attention to scaphoid: 1mm small foc scaphoid reformats, parallel to the long axis of the scaphoid bone.
<b>Comments:</b>	Recon 1 is a single thin helical group of the wrist/hand for direct mpr. Recon 2 is the wrist/hand, bone algorithm ct going to PACS. Recon 3 is wrist/hand, standard algorithm ct going to PACS.
<b>Images required in PACS</b>	Scouts, 1.25mm x 1.25mm axial wrist/hand bone, 1.25mm x 1.25mm axial wrist/hand standard, 2mm x 2mm sagittal wrist/hand, 2mm x 2mm coronal wrist/hand, Dose Report