

**RIH - SHOULDER CT
SIEMENS DEFINITION AS+ PROTOCOL**

Indication: fracture, dislocation, osteomyelitis, bone injury, bone tumor.

Position/Landmark	Supine , feet first Zero Appropriately																														
Topogram Direction	Craniocaudal																														
Respiratory Phase	Suspension																														
Scan Type	Helical																														
Ref kV/Ref mAs/Rotation time (sec) Pitch / Speed (mm/rotation) Safire Strength / Dose Optimization	Care kV 120 / Care Dose4D 100 / 1 sec .8:1 , 32.00mm 3 / 4																														
Detector width x Rows = Beam Collimation	0.625mm x 64 = 40mm (128 x .6mm)																														
Average Tube Output	ctdi – 11.0mGy dlp – 280mGy.cm																														
Helical Set Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>axial soft shoulder</td> <td>3mm x 3mm</td> <td>I40s medium</td> <td>pac</td> </tr> <tr> <td>2</td> <td>axial bony shoulder</td> <td>3mm x 3mm</td> <td>I70h very sharp</td> <td>pac</td> </tr> <tr> <td>3</td> <td>coronal shoulder</td> <td>3mm x 3mm</td> <td>I70h very sharp</td> <td>pac</td> </tr> <tr> <td>4</td> <td>sagittal shoulder</td> <td>3mm x 3mm</td> <td>I70h very sharp</td> <td>pac</td> </tr> <tr> <td>5</td> <td>thin shoulder</td> <td>.75mm x .7mm</td> <td>I70h very sharp</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination	1	axial soft shoulder	3mm x 3mm	I40s medium	pac	2	axial bony shoulder	3mm x 3mm	I70h very sharp	pac	3	coronal shoulder	3mm x 3mm	I70h very sharp	pac	4	sagittal shoulder	3mm x 3mm	I70h very sharp	pac	5	thin shoulder	.75mm x .7mm	I70h very sharp	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination																											
1	axial soft shoulder	3mm x 3mm	I40s medium	pac																											
2	axial bony shoulder	3mm x 3mm	I70h very sharp	pac																											
3	coronal shoulder	3mm x 3mm	I70h very sharp	pac																											
4	sagittal shoulder	3mm x 3mm	I70h very sharp	pac																											
5	thin shoulder	.75mm x .7mm	I70h very sharp	terarecon																											
Scan Start / End Locations	determined by technologist or radiologist to include the anatomy of interest																														
DFOV	18cm decrease appropriately																														
IV Contrast Volume / Type / Rate	75mL Iohexol (Omnipaque 350) / 2mL per second if needed																														
Scan Delay	65 seconds																														
2D/3D Technique Used	Workstream 4D mpr of 3mm x 3mm coronal and sagittal shoulder series (auto-batch off), average mode, auto-transferred to PACS																														
Comments: Recon 5 is a thin helical volume of the shoulder that is archived to the TeraRecon server.																															
Images required in PACS	Topograms, 3mm x 3mm axial shoulder bone, 3mm x 3mm axial shoulder standard, 3mm x 3mm sagittal shoulder, 3mm x 3mm coronal shoulder, Patient Protocol																														