

## RIH - ELBOW CT GE LIGHTSPEED VCT PROTOCOL

**Indication: fracture, dislocation, osteomyelitis, bone injury, bone tumor.**

<b>Position/Landmark</b>	Supine , feet first Zero Appropriately			
<b>Topogram Direction</b>	Craniocaudal			
<b>Respiratory Phase</b>	Any			
<b>Scan Type</b>	Helical			
<b>KV / mA / Rotation time (sec)</b> <b>Pitch / Speed (mm/rotation)</b> <b>Noise Index / ASiR / Dose Reduction</b>	120kv / smart mA (100-450) / 0.5 sec 0.984:1 , 39.37mm 16.0 / 20 / 20%			
<b>Detector width x Rows = Beam Collimation</b>	0.625mm x 64 = 40mm			
<b>Helical Set</b>	recon	body part	thickness/ spacing	recon destination .
Slice Thickness/ Spacing			algorithm	
Algorithm	1	<b>elbow bone</b>	2.5mm x 2.5mm	bone pacs
Recon Destination	2	thin elbow	.6mm x .6mm	bone for dmpr
	3	<b>elbow soft tissue</b>	2.5mm x 2.5mm	standard pacs
<b>Scan Start / End Locations</b>	determined by technologist or radiologist to include the anatomy of interest			
<b>DFOV</b>	18cm decrease appropriately			
<b>IV Contrast Volume / Type / Rate</b>	75mL Iohexol (Omnipaque 350) / 2mL per second if needed			
<b>Scan Delay</b>	65 seconds			
<b>2D/3D Technique Used</b>	DMPR of 3mm x 3mm coronal and sagittal elbow series (auto-batch off), average mode, auto-transferred to PACS  Also, there is a 3mm x 3mm true axial reformat if needed due to the patient's position.			
<b>Comments:</b>	Recon 1 is the 2.5mm x 2.5mm elbow, bone algorithm ct going to PACS. Recon 2 is a single thin helical group of the elbow for direct mpr. Recon 3 is the 2.5mm x 2.5mm elbow, standard algorithm ct going to PACS.			
<b>Images required in PACS</b>	Scouts, 2.5mm x 2.5mm axial elbow bone, 2.5mm x 2.5mm axial elbow standard, 3mm x 3mm sagittal elbow, 3mm x 3mm coronal elbow, Dose Report			