

**RIH – TOTAL BODY ANGIO / SOLID ORGAN INJURY FOR TRAUMA
SIEMENS DEFINITION AS20 PROTOCOL**

Position/Landmark	Head first or feet first-Supine 1cm superior to skull vertex																																								
Topogram Direction	Craniocaudal / Craniocaudal																																								
Scan Type	Helical																																								
Ref kV/Ref mAs/Rotation time (sec) Pitch / Speed (mm/rotation) Safire Strength / Dose Optimization	5sec Care kV 120 / Care Dose4D 210 / 0.5sec 1.2:1 , 15.00mm 3 / 6																																								
Detector width x Rows = Beam Collimation	0.625mm x 20 = 12.5mm																																								
Average Tube Output	11.0mGy dlp – 1456 mGy.cm																																								
Helical Set Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1"> <thead> <tr> <th></th> <th>body part</th> <th>thickness/ spacing</th> <th>recon algorithm</th> <th>recon destination</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>total body ct angio</td> <td>2mm x 2mm</td> <td>B30s medium</td> <td>pac</td> </tr> <tr> <td>2</td> <td>axial chest abd pelvis</td> <td>5mm x 5mm</td> <td>I40f medium</td> <td>pac</td> </tr> <tr> <td>3</td> <td>axial t and l spine</td> <td>3mm x 3mm</td> <td>I70f very sharp</td> <td>pac</td> </tr> <tr> <td>4</td> <td>lungs</td> <td>5mm x 5mm</td> <td>I70f very sharp</td> <td>pac</td> </tr> <tr> <td>5</td> <td>coronal chest abd pelvis</td> <td>5mm x 5mm</td> <td>I40f medium</td> <td>pac</td> </tr> <tr> <td>6</td> <td>thin body ct angio</td> <td>.75mm x .7mm</td> <td>B30s medium</td> <td>mpr</td> </tr> <tr> <td>7</td> <td>thin t and l spine</td> <td>.75mm x .7mm</td> <td>I70f very sharp</td> <td>mpr</td> </tr> </tbody> </table>		body part	thickness/ spacing	recon algorithm	recon destination	1	total body ct angio	2mm x 2mm	B30s medium	pac	2	axial chest abd pelvis	5mm x 5mm	I40f medium	pac	3	axial t and l spine	3mm x 3mm	I70f very sharp	pac	4	lungs	5mm x 5mm	I70f very sharp	pac	5	coronal chest abd pelvis	5mm x 5mm	I40f medium	pac	6	thin body ct angio	.75mm x .7mm	B30s medium	mpr	7	thin t and l spine	.75mm x .7mm	I70f very sharp	mpr
	body part	thickness/ spacing	recon algorithm	recon destination																																					
1	total body ct angio	2mm x 2mm	B30s medium	pac																																					
2	axial chest abd pelvis	5mm x 5mm	I40f medium	pac																																					
3	axial t and l spine	3mm x 3mm	I70f very sharp	pac																																					
4	lungs	5mm x 5mm	I70f very sharp	pac																																					
5	coronal chest abd pelvis	5mm x 5mm	I40f medium	pac																																					
6	thin body ct angio	.75mm x .7mm	B30s medium	mpr																																					
7	thin t and l spine	.75mm x .7mm	I70f very sharp	mpr																																					
Scan Start / End Locations DFOV	1cm superior to skull through the toes 38cm decrease appropriately																																								
IV Contrast Volume / Type / Rate	Hand inject 30mL Iohexol (Omnipaque 350) then wait 5 minutes 80mL Iohexol (Omnipaque 350) / 3mL per second 15mL saline / 3mL per second start monitoring scans at this point 100mL Iohexol (Omnipaque 350) / 4mL per second 40mL saline / 4mL per second																																								
Scan Delay	bolus tracking at aortic arch, the trigger is +150 HU																																								
2D/3D Technique Used	2mm x 2mm left sagittal/oblique carotid, 2mm x 2mm right sagittal/oblique carotid, 2mm x 2mm coronal carotids, 2mm x 2mm axial small fov carotids, 5mm x 5mm coronal chest abdomen pelvis, 3mm x 3mm coronal lower extremities. Recon 7 is a thin bone volume incase spines are ordered.																																								
Comments: 80cc of contrast is pre injected to be able to evaluate for solid organ injury. The second 100cc is injected immediately after the first injection is completed to perform a smart prepped total body ct angiogram.																																									
<ol style="list-style-type: none"> 1. Hand inject 30mL Iohexol (Omnipaque 350) 2. Perform scout images 3. Plan all three recons 4. Perform pre-monitoring scan at the aortic arch, Place the ROI in the arch 5. Inject the 80mL Iohexol (Omnipaque 350) pre-scan contrast 																																									
As soon as the power injector hold begins, inject the CTA contrast and begin the monitoring phase ct angiogram.																																									
Images required in PACS	Topograms, 2mm x 2mm axial cta chest to toes, 2mm x 2mm axial carotid cta, 2mm x 2mm left sagittal/oblique carotid, 2mm x 2mm right sagittal/oblique carotid, 2mm x 2mm coronal carotids, 5mm x 5mm coronal chest abdomen pelvis, 5mm x 5mm axial lungs, 3mm x 3mm axial coronal and sagittal cervical, thoracic and lumbar spines, 3mm x 3mm coronal lower extremity cta, Patient Protocol																																								